

Digital Dental Labs - Scanner Connection Guide

3Shape TRIOS® 1. Visit us.3shapecommunicate.com 2. Search digitaldentalportal@gmail.com 3. Select Digital Dental Labs	Medit 1. Register 2. Add Partner 3. Upload STL/PLY/OBJ
iTero™ 1. Find a Laboratory 2. Company ID: 506826 3. Select Digital Dental Labs	Planmeca® 1. Find a Lab 2. Search 3. Add Lab 4. Submit
Carestream 1. Email the lab 2. Search Digital Dental Labs 3. Submit STL files	Shining 3D 1. Dental Cloud 2. Connections 3. Confirm 4. Send scans
CEREC® 1. Log in 2. My CEREC Connect 3. Search/Add Lab 4. Submit case	Dexis 1. ScanFlow 2. DentalConnect 3. Add Partner 4. Request connection